



marLED[®] E3

Top-quality examination lights

**100,000
Lux**



marLED® E3 Top-quality examination lights

State-of-the-art LED technology will guarantee highest light output without getting you hot under the collar. Thanks to the easy to move swivel mounting of the lamp, reliable positioning is easy wherever you need light.

Excellent color reproduction, color temperature designed to mimic natural daylight and a long service life combined with low operating costs make this examination light the perfect tool for your day-to-day work.

Ideal for clinics and surgeries

- in outpatient departments
- in prep and recovery rooms
- in dental and oral surgery
- in oral and maxillofacial surgery
- in veterinary medicine
- in day clinics and outpatient clinics


Performance. Efficiency. Technology.
Light for the most demanding requirements



Operation –
is brilliantly simple!



Safe hygiene –
with no corners or edges



TIR optics for brilliant light –
accurately calculated



Reliable positioning –
lightweight and easy to move



Technology at a glance

	marLED® E3
Electrical supply voltage provided via integrated wide-range power supply	100 – 240 VAC/50-60 Hz
Color temperature	4,300 K
Illuminance at a distance of 1.0 m	100,000 Lux
Light field diameter	170 mm
Power consumption of lights	3 x 5 watt
Anticipated service life	> 60,000 h

Mounting variants



Mobile version



Ceiling version



Wall version

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co. Ltd.

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 505 7838
malaysia@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

London · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
P.O. Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

